

Citrate Synthase (CS) Activity Assay Kit

Note: Take two or three different samples for prediction before test.

Detection equipment: Spectrophotometer

Cat No:NA0792

Size:50T/24S

Components:

Extract solution:25 mL×1, store at -20°C.

Reagent I:5 mL×1, store at 4°C.

Reagent II:0.3 mL×1, store at -20°C.

Reagent III: 90 mL×1, store at 4°C.

Reagent IV: 4 mL×1, store at 4°C.

Reagent V: Powder×4, store -20°C. Add 500 µL of distilled water when the solution will be used. It is suggested that the unused reagent should still be storage at -20°C.

Reagent VI: Powder×2, store at -20°C. Add 1.5 mL of distilled water when the solution will be used. It is suggested that the unused reagent should still be storage at -20°C.

Description:

Citrate Synthase (CS, EC 2.3.3.1) is widely exists in animals, plants, microorganism and mitochondrial matrix of cultured cells. It is the first rate-limiting enzyme in the tricarboxylic acid cycle and one of the main regulatory sites.

CS catalyzes acetyl CoA and acetoacetic acid to generate citryl coenzyme A, further hydrolysis to produce citric acid. The reaction promoted the transformation of colorless DTNB to yellow TNB, which has absorption at 412 nm.

Required but not provided:

Spectrophotometer, low temperature centrifuge, water bath, mortar/homogenizer, adjustable pipette, 1 mL glass cuvette and distilled water.

Protocol:

I. Sample Extraction:

Isolation of cytoplasmic proteins and mitochondrial proteins from tissues, bacteria, or cells:

1. Take 0.1 g of tissue or 5 million cells, add 1 mL of Extract solution and 10 µL of Reagent II, homogenate with homogenizer.
2. Centrifuge at 600 ×g and 4°C for 5 minutes.
3. Take supernatant to another centrifuge tube, centrifuge at 11000 ×g and 4°C for 10 minutes.
1. 4. The supernatant is a plasma extract that can be used to determine CS leakage from mitochondria.
5. Add 200 µL of Reagent I and 2 µL of Reagent II into precipitate, mix thoroughly to detect the activity of

CS and the detection of protein concentration.

II. Procedure

1. Preheat spectrophotometer for 30 minutes, adjust wavelength to 412 nm, set zero with distilled water.
2. Preheat Reagent III at 25°C (general species) or 37°C (mammals) water bath for 10 minutes (guaranteed no precipitation).
3. Procedure test

Reagent name (μL)	Test tube (T)	Control tube (C)
Reagent III	860	930
Reagent IV	35	35
Reagent V	35	-
Sample	35	35
Reagent VI	35	-

Add reagents orderly to the 1 mL glass cuvette, record the time when adding Reagent VI, record the absorbance A1 of 10s at 412 nm. Then place the cuvette with reaction solution to 37°C or 25°C water bath for 2 min. Take out and wipe to dry the cuvette, record the absorbance A2 of 412 nm at 130s, test tube and control tube all need detect.

Test tube: $\Delta A_1 = A_2 - A_1$, Control tube $\Delta A_1' = A_2 - A_1$, $\Delta A = \Delta A_1 - \Delta A_1'$.

III. Calculation

1. Tissue

(1). Protein concentration

Unit definition: One unit of enzyme activity is defined as the amount of enzyme catalyzes the generation of 1 nmol TNB at 37°C or 25°C every milligram of tissue protein per minute.

$$CS(U/mg \text{ prot}) = \Delta A \div (\epsilon \times d) \times V_{RT} \div (C_{pr} \times V_s) \div T = 1050 \times \Delta A \div C_{pr}$$

ϵ : Molar extinction coefficient of TNB, $13.6 \times 10^{-3} \text{ mL}/(\text{nmol} \cdot \text{cm})$;

V_{RT} : Reaction volume, 1 mL;

d : Cuvette diameter (cm), 1 cm;

V_s : Sample volume, 0.035 mL;

T : Reaction time (min), 2 minutes;

C_{pr} : Protein concentration after precipitation dissolution, mg/mL.

Note:

1. Samples and all reagents place on ice, in order to avoiding denaturation and lose activity.
2. The reaction solution of cuvette should place 37°C or 25°C. Add a certain amount of distilled water to a small beaker, then the small beaker place in water bath at 37°C or 25°C. Place the cuvette with the solution in the beaker during the reaction.
3. Two people do this experiment at the same time, one person colorimetric, the other person timing to ensure the accuracy of experiment results.
4. It is recommended to use the concentration of sample protein to calculate the enzyme activity. If the fresh weight of sample is used to calculate, the enzyme activity of enzyme solution extracted from

cytoplasm should be measured. The sum of enzyme activity in supernatant and precipitation is the total enzyme activity.

5. Appendix: Calculation formula of sample weight.

Unit definition: One unit of enzyme activity is defined as the amount of enzyme catalyzes the generation of 1 nmol TNB at 37°C or 25°C every gram tissue per minute in the reaction system.

$$CS_S(\text{U/mg prot}) = \Delta A_S \div \epsilon \div d \times V_{RT} \div (W \times V_S \div V_E) \div T = 1061 \times \Delta A_S \div W$$

$$CS_P(\text{U/mg prot}) = \Delta A_P \div \epsilon \div d \times V_{RT} \div (W \times V_S \div V_P) \div T = 212 \times \Delta A_P \div W$$

$$CS(\text{U/mg prot}) = CS_S + CS_P = 1061 \times \Delta A_S \div W + 212 \times \Delta A_P \div W$$

ΔA_1 : The measured value of supernatant;

ΔA_2 : The measured value of precipitation;

ϵ : Molar extinction coefficient of TNB, $13.6 \times 10^3 \text{ mL}/(\text{nmol} \cdot \text{cm})$;

V_{RT} : Reaction volume, 1 mL;

d : Cuvette diameter(cm), 1 cm;

V_S : Sample volume, 0.035 mL;

V_E : Extract solution volume, 1.01 mL;

V_P : Total volume of precipitation, 0.202 mL;

T : Reaction time(min), 2 minutes.

W : Sample weight, g.

Experimental instances:

1. Add 1mL of Extract solution and 10 μ L of Reagent II into 0.1g of mouse heart, homogenate and grind. Take supernatant and centrifuge, take the supernatant and sediment. Add 200 μ L of Reagent I and 2 μ L of Reagent II to the sediment, test according to the measured steps. Calculate in the supernatant: $\Delta A_1 = A_2 - A_1 = 1.213 - 0.7 = 0.513$, $\Delta A_1' = A_2' - A_1' = 0$, $\Delta A_S = \Delta A_1 - \Delta A_1' = 0.513$, Calculate in the precipitation: $\Delta A_1 = A_2 - A_1 = 0.447 - 0.133 = 0.314$, $\Delta A_1' = A_2' - A_1' = 0$, $\Delta A_P = \Delta A_1 - \Delta A_1' = 0.314$, calculate the enzyme activity according to sample weight: $CS(\text{U/g weight}) = CS_S + CS_P = 1061 \times \Delta A_S \div W + 212 \times \Delta A_P \div W = 1061 \times 0.513 \div 0.1 + 212 \times 0.314 \div 0.1 = 6108.7 \text{ U/g weight}$.

Recent Product citations

[1] Ming Song, Fangfang Chen, Yihui Li, et al. Trimetazidine restores the positive adaptation to exercise training by mitigating statin-induced skeletal muscle injury. *Journal of Cachexia, Sarcopenia and Muscle*. November 2017; (IF10.754)

[2] Zhang J, Lv J, Xie J, et al. Nitrogen Source Affects the Composition of Metabolites in Pepper (*Capsicum annuum* L.) and Regulates the Synthesis of Capsaicinoids through the GOGAT-GS Pathway[J]. *Foods*, 2020, 9(2): 150.

Reference:

[1] Agostinho F R, Réus G Z, Stringari R B, et al. Treatment with olanzapine, fluoxetine and olanzapine/fluoxetine alters citrate synthase activity in rat brain[J]. *Neuroscience letters*, 2011, 487(3):

278-281.

Related products:

NA0843/NA0601	CoenzymeINAD(H) Content Assay Kit
NA0795/NA0554	NAD Kinase (NADK) Assay Kit
NA0818/NA0576	NADH oxidase (NOX) Activity Assay Kit
NA0785/NA0544	NAD Malic Enzyme (NAD-ME) Assay Kit