

## 5'- Nucleotidase (5'-NT) Activity Assay Kit

**Note:** It is necessary to predict 2-3 large difference samples before the formal determination.

**Operation Equipment:** Spectrophotometer/ Microplate reader

**Cat No:** NA0262

**Size:** 50T/24S

### Components:

Extracting solution: Liquid 30 mL×1. Storage at -20°C.

Reagent I: Powder ×2. Storage at -20°C.

Reagent II: Liquid 12 mL×2. Storage at 4°C.

Reagent III: Liquid 30 mL×1. Storage at 4°C.

Reagent IV: Liquid 25 mL×1. Storage at 4°C.

Reagent V: Powder ×1. Storage at 4°C. Before use, add 12 mL of distilled water, fully dissolve, and store the unused reagent at 4°C for two weeks.

Reagent VI: Powder×1. Storage at 4°C. Before use, add 12 mL of distilled water, fully dissolve, and store the unused reagent at 4°C for two weeks.

Reagent VII: Liquid 12 mL×1. Storage at room temperature.

Standard solution: Powder×1. Storage at 4°C. 8 mg of phosphorus standard. Before use, 4.6 mL of Reagent IV is added to prepare a standard solution of 10 μmol/mL. After dissolution, the solution is stored at 4°C.

Working solution: Reagent I are added into a bottle of Reagent II to dissolve completely; the unused reagents are packed and stored at - 20°C for one week, and prepare when the solution will be used.

Preparation of phosphorus determination reagent: prepare according to the proportion of H<sub>2</sub>O: Reagent V: Reagent VI: Reagent VII = 2:1:1:1, and the prepared phosphorus determination reagent shall be light yellow. If colorless, reagent fails; if blue, it is phosphorus pollution (please use how much to match as required).

### Product Description:

5'-nucleotidase (5'-NT) is a kind of hydrolase with low substrate specificity, which can act on a variety of nucleotides. It widely exists in various plant, animal tissues, serum and plasma. 5'-NT is a special phosphate hydrolase, which acts on nucleoside-5'-phosphate such as AMP (adenosine-5'-phosphate or adenosine monophosphate) to produce inorganic phosphate and nucleoside. The activity of 5'-NT can be calculated by determining the content of inorganic phosphorus.

### Reagents and Equipment Required but Not Provided:

Balance, Spectrophotometer, desktop centrifuge, cryogenic centrifuge, constant temperature water bath/constant temperature incubator, 1 mL glass cuvette, transferpettor, mortar/homogenizer, ice, distilled water.

## Procedure:

**I. Sample preparation** (the sample size can be adjusted appropriately, and the specific proportion can be referred to the literature):

1. Tissue: The ratio of mass (g): volume of Extracting solution (mL) is 1:5-10 (it is recommended to weigh about 0.1 g and add 1 mL of Extracting solution), homogenize on ice, centrifuge at 4°C, 15000 g for 10 min, and place the supernatant on ice for testing.

2. Cells: The ratio of the number of cells ( $10^4$ ): the volume of distilled water (mL) is 500-1000:1 (it is recommended to add 1 mL distilled water to 5 million cells), the cells are broken by ice bath ultrasonic wave (power 300W, ultrasonic 3s, interval 7s, total time 3 min); then the cells are centrifuged at 4°C, 15000g for 10 min, and the supernatant is put on ice for testing.

3. Liquid: direct detection.

## II. Determination procedure:

1. Preheat the Spectrophotometer for 30 minutes, adjust the wavelength to 660 nm, set zero with distilled water.
2. The starch standard solution is diluted with Reagent IV to 0.48、0.24、0.12、0.06、0.03、0.015  $\mu\text{mol/mL}$ .
3. Add reagents with the following list: (Operate in 1.5 mL EP tube)

### (1) Enzymatic reaction

| Reagent ( $\mu\text{L}$ )                                                                        | Test tube | Control tube |
|--------------------------------------------------------------------------------------------------|-----------|--------------|
| Sample                                                                                           | 100       | 100          |
| Working solution                                                                                 | 400       |              |
| ortex mixing, 37°C (mammalian) or 25°C (plant and other) reaction for 30 min                     |           |              |
| Reagent III                                                                                      | 500       | 500          |
| Working solution                                                                                 | -         | 400          |
| Vortex mixing, 25°C, 8000 rpm centrifugation for 10 min, take the supernatant for color reaction |           |              |

### (2) Color reaction

| Reagent ( $\mu\text{L}$ )        | Test tube | Control tube | Standard tube | Blank tube |
|----------------------------------|-----------|--------------|---------------|------------|
| Supernatant                      | 400       | 400          | -             | -          |
| Standard                         | -         | -            | 400           | -          |
| Reagent IV                       | -         | -            | -             | 400        |
| Phosphorus determination reagent | 800       | 800          | 800           | 800        |

Vortex mixing, 40°C color for 10 min; take 1 mL of reaction solution in 1 mL glass cuvette, measure the absorbance value A at 660 nm, respectively record as  $A_T$ ,  $A_C$ ,  $A_S$ ,  $A_B$ , calculate  $\Delta A_S = A_S - A_B$ ,  $\Delta A_T = A_T - A_C$  (blank tube only needs to measure 1-2 times).

## III. Calculation:

1. Drawing of standard curve: draw the standard curve with  $\Delta A_S$  as y axis, and the standard solution concentration as x axis, and get the standard equation  $y=kx+b$ , and bring the  $\Delta A$  into the equation to get

x( $\mu\text{mol/mL}$ ).

## 2. Calculation of 5'-NT activity

### (1) Calculated according to protein concentration

Unit definition: One unit of enzyme is defined as the amount of enzyme catalyzes the production of 1 nmol inorganic phosphorus per minute every milligram tissue protein in the reaction system.

$$5'\text{-NT activity (U/mg prot)} = x \times V_{RT} \div (V_S \times C_{pr}) \div T \times 10^3 = 333.3 \times x \div C_{pr}$$

### (2) Calculated by sample mass

Unit definition: One unit of enzyme is defined as the amount of enzyme catalyzes the production of 1 nmol inorganic phosphorus per minute every milligram tissue in the reaction system.

$$5'\text{-NT activity (U/g mass)} = x \times V_{RT} \div (W \times V_S \div V_{ST}) \div T \times 10^3 = 333.3 \times x \div W$$

### (3) Calculated by cell number

Unit definition: One unit of enzyme is defined as the amount of enzyme catalyzes the production of 1 nmol inorganic phosphorus per minute every  $10^4$  cells in the reaction system.

$$5'\text{-NT activity (U}/10^4 \text{ cell)} = x \times V_{RT} \div (\text{cell number} \times V_S \div V_{ST}) \div T \times 10^3 = 333.3 \times x \div \text{cell number}$$

### (4) Calculated according to volume of liquid:

Unit definition: One unit of enzyme is defined as the amount of enzyme catalyzes the production of 1 nmol inorganic phosphorus per minute every milliliter liquid in the reaction system.

$$5'\text{-NT activity (U/mL)} = x \times V_{RT} \div V_S \div T \times 10^3 = 333.3 \times x$$

$V_S$ : sample volume added in enzymatic reaction, 0.1 mL;  $V_{RT}$ : total volume of enzymatic reaction, 1 mL;  $V_{ST}$ : volume added in Extracting solution, 1 mL; W: sample mass, g; Cpr: sample protein concentration, mg/mL; cell number: in tens of thousands; T: enzymatic reaction time, 30 min;  $10^3$ : unit conversion, 1  $\mu\text{mol} = 10^3 \text{ nmol}$ .

## Note:

When the absorbance value is greater than 1 or  $\Delta A$  is greater than 1, it is suggested that the sample be diluted with Reagent IV before determination.

## Experimental example:

1. Take 0.1 g of mouse liver, and then take the sample for treatment. take the supernatant and operate according to the determination steps. Calculate:  $\Delta A_T = A_T - A_C = 0.723 - 0.534 = 0.189$ , and bring the standard curve  $y = 2.3928x + 0.0165$ , calculate  $x = 0.0721$ , calculate the enzyme activity according to the sample quality:

$$5'\text{-NT activity (U/g mass)} = 333.3 \times x \div W = 333.3 \times 0.0721 \div 0.1 = 240.31 \text{ U/g mass.}$$

2. Take 0.1 g of barnyard grass for sample treatment. take the supernatant and operate according to the determination steps. Calculate:  $\Delta A_T = A_T - A_C = 0.367 - 0.281 = 0.086$ , and bring in the standard curve  $y = 2.3928x + 0.0165$ , calculate  $x = 0.0290$ , calculate the enzyme activity according to the sample quality:

$$5'\text{-NT activity (U/g mass)} = 333.3 \times x \div W = 333.3 \times 0.0290 \div 0.1 = 96.657 \text{ U/g mass.}$$

## Related products:

NA0784/NA0543 Creatine Kinase (CK) Activity Assay Kit

NA0288/NA0287 Pyrroline-5-carboxylic Acid Synthase (P5CS) Activity Assay Kit

NA0738/NA0495 Laccase Activity Assay Kit

NA0724/NA0482 Isocitrate Lyase (ICL) Activity Assay Kit

NA0640/NA0377 Acetate Kinase (ACK) Activity Assay Kit